


**Lab 2**

## Statistical Copolymers via Free-Radical Copolymerization

**Report Due:** *In Lecture*, Monday, February 25**Report Revision Due:** *In Lecture*, Monday, March 11**Introduction**

Copolymerization of multiple monomers allows for the synthesis of new polymeric materials with properties that are often substantially different from homopolymers made from the individual monomers. This experiment will allow you to perform the radical copolymerization of styrene and methylmethacrylate and determine the relative reactivities of these monomers in this process. As in Lab 1, polymerizations will be initiated via the thermal decomposition of benzoyl peroxide to benzoyl radicals. You will also learn to characterize polymers by nuclear magnetic resonance (NMR) spectroscopy.

As early as the 1930's it was recognized that in the polymerization of equimolar mixtures of two different monomers, the tendencies of the monomers to incorporate into the growing chain will differ. The result of the variable incorporation of monomer is that the composition of the feed (and thus the copolymer) will change as the polymerization proceeds. Variable monomer incorporation arises from selectivity of monomer addition at the growing chain end. The copolymerization of monomers  $M_1$  and  $M_2$  can give rise to four different propagation opportunities (Figure 2-1):

**Figure 2-1.**

Here,  $k_{11}$  and  $k_{22}$  are the rate constants for self-propagation, and  $k_{12}$  and  $k_{21}$  are the rate constants for cross-propagation. The relative reactivity ratios for  $M_1$  and  $M_2$  in this polymerization can be described as  $r_1 = (k_{11}/k_{12})$  and  $r_2 = (k_{22}/k_{21})$ , respectively.

The rate of consumption of  $M_1$  initially can be described by:

$$(2.1) \quad -\frac{d[M_1]}{dt} = k_{11}[M_1][M_1\bullet] + k_{21}[M_1][M_2\bullet]$$

And of  $M_2$ :

$$(2.2) \quad -\frac{d[M_2]}{dt} = k_{22}[M_2][M_2\bullet] + k_{12}[M_2][M_1\bullet]$$

The consumption of  $M_1$  and  $M_2$  is due to their incorporation into the growing copolymer. By writing the equation above as a ratio of rates for the two monomers, and by applying the steady state approximation, a ratio of monomer incorporation into the copolymer can be expressed as the instantaneous copolymerization equation:

$$(2.3) \quad \frac{d[M_1]}{d[M_2]} = \frac{[M_1](r_1[M_1] + [M_2])}{[M_2]([M_1] + r_2[M_2])}$$

If the monomer feed concentrations of  $M_1$  and  $M_2$  (mole fractions) are denoted  $f_1$  and  $f_2$ , and the mole fraction of these monomers that are incorporated into the copolymer as  $F_1$  and  $F_2$ , the instantaneous copolymerization equation can be defined for  $F_1$ :

$$(2.4) \quad F_1 = \frac{r_1 f_1^2 + f_1 f_2}{r_1 f_1^2 + 2 f_1 f_2 + r_2 f_2^2}$$

Experimentally, reactivity ratios are often determined by performing a series of polymerizations under identical experimental conditions, but beginning each polymerization at different concentrations of the two monomers. The polymerizations are terminated at low conversion (3-10%), and the resulting polymer is characterized to determine the levels of incorporation of the two monomers ( $F_1$  and  $F_2$ ).

A variety of methods have been devised for extracting  $r_1$  and  $r_2$  from copolymer feed ( $f_n$ ) and composition ( $F_n$ ) data. In the days before computerized data analysis, Fineman and Ross<sup>1</sup> and Kelen and Tüdös<sup>2</sup> developed forced linearization methods for dealing with the complexity of equation 2.4. Today, non-linear least-squares curve-fitting programs (such as Kaleidagraph) allow for the determination of  $r_1$  and  $r_2$  without resorting to a forced linearization protocol.<sup>3</sup> Values for  $r_1$  and  $r_2$  can be obtained by a form of the instantaneous copolymerization equation given in terms of  $f_1$  and  $F_1$  by substituting  $(1 - f_1)$  for  $f_2$  in equation 2.4:

$$(2.5) \quad F_1 = \frac{r_1 f_1^2 + f_1(1 - f_1)}{r_1 f_1^2 + 2 f_1(1 - f_1) + r_2(1 - f_1)^2}$$

<sup>1</sup> Fineman, M.; Ross, S.D. *J. Polym. Sci.*, **1950**, 5, 259.

<sup>2</sup> Kelen, T.; Tüdös, F. *J. Macromol. Sci.-Chem.*, **1975**, A9, 1.

<sup>3</sup> For a more detailed discussion, including references on how to implement this fit, see: (a) Bauduin, G.; Boutevin, B.; Belbachir, M.; Meghabar, R. *Macromolecules* **1995**, 28, 1750, and references therein. (b) G. Odian, *Principles of Polymerization*, pp. 467-470.

In this case, the data of a plot of  $f_1$  vs.  $F_1$  can be analytically fit to equation 2.5 to obtain values for  $r_1$  and  $r_2$ .

## **Experimental**

### **Bulk copolymerization of styrene and *n*-butylmethacrylate (Feb 5/7)**

*The goal of this experiment is to observe the consequence of copolymer ratio on the physical properties of a copolymer material.*

- In a 6" test tube, pour a total of 5 mL of styrene and *n*-butylmethacrylate; feel free to choose any ratio of the two monomers, though styrene content between 5 and 25 wt% will give you a nice spectrum of plastic flexibility. Add 0.05 g of benzoyl peroxide, and swirl to dissolve. Follow the instructions in Lab 1 for spraying an aluminum pan with silicone mold release and polymerizing the mixture in the 70 °C oven.
- *Next week:* Remove your polymerized sample from the aluminum pan. How does this sample compare with the homopolymers from Lab 1?

### **Copolymerization of styrene and methylmethacrylate (Feb 5/7)**

*The goal of this experiment is to copolymerize styrene and methyl methacrylate using a free radical initiator, and determine reactivity ratios for the resulting copolymers.*

- In three separate 20-mL scintillation vials, put 5 mL (total) of a styrene-methyl methacrylate monomer mixture. For best accuracy, we recommend that you use balances in the hoods to measure reagent masses rather than working with reagent volumes. Use the ratios of styrene to methylmethacrylate that were listed in your pre-lab handout. Add 0.1 g of benzoyl peroxide to each of the tubes and swirl to dissolve. Repeat the procedure outlined for the styrene polymerization performed in Lab 1. Make sure to record the polymerization time and temperature, and dry your product in the vacuum oven as you did with the styrene homopolymer.

### **NMR analysis of copolymer ratios (Feb 12/14)**

- Determine the mass yield of your samples. If the yield is greater than 10%, or you isolate less than 50 mg of polymer (about 1%), you should disregard the data from that polymerization run.
- Check out three NMR tubes from the stockroom. Put approximately 25 mg of each of your polymer samples into a small vial along with about 1 mL of CDCl<sub>3</sub> to dissolve, and allow to sit for about 30 minutes to fully dissolve. Then, transfer each solution to an NMR tube. CDCl<sub>3</sub> is very expensive, so don't spill or waste it. Put your samples in the class NMR tube rack according to the TA's instructions, and label each tube with your group number and the

fraction of styrene ( $f_{\text{styrene}}$ ) monomer you used. (For those of you familiar with taking NMR spectra of organic samples, 25 mg will seem like too much material. However, because peaks in NMR spectra of polymers tend to be broad, you'll need more than you're used to.)

Your NMR spectra for this class will be acquired by the TA's in the Chemistry Department NMR Facility, but you will be expected to work up the data yourself. Raw FID's will be available for processing in the data directories of the Polymer Lab NMR account. Each FID file will be titled according to the date and time it was acquired; a listing of which files are yours will be posted on the "Data" page of the course website.

You'll be working up your NMR data using ACDLabs' NMR Processor for PC. Instructions on how to use this program are in Appendix A.

You should use the integrations you performed on the NMR spectra to determine  $F_1$  and  $F_2$  for your copolymer samples. Remember to factor in the number of protons that contribute to each NMR peak(s) per molecule. As part of an Assignment, you will be expected to post your data to the course WebCT Discussion Board.

### **Lab Report** (due in lecture, Monday, February 25)

Prepare a lab report titled, "Determination of reactivity ratios in the free radical polymerization of styrene and methylmethacrylate." Using nonlinear least-squares fitting to Equation 2.5, calculate  $r_1$  and  $r_2$  from the class'  $f_1$  and  $F_1$  data for the copolymerization of styrene and methylmethacrylate. You can use Kaleidagraph, a program available in the Microcomputer Lab, to do this; or, you can use any program capable of non-linear least-squares fitting. Instructions on how to use Kaleidagraph are in Appendix B at the end of this handout.

The report should be single-spaced, and any figures, tables or graphs should be inserted in the text. Your report should be structured like a journal article in an American Chemical Society journal (such as *Macromolecules*), like the one you examined as part of Assignment 4. As you did in that assignment, you may want to consult websites or articles on how to write a scientific paper.

Please be concise in your writing. Be sure to list the names of all of your group members as coauthors; as the corresponding author, your name should be designated with a \*. At any point, if you wonder how to structure or phrase something, look at a couple of *Macromolecules* papers and copy their format.

Overall, your report should stand on its own, independent of this course. If a reader from another university with no knowledge of your assignment were to read your report, would they understand it? If not, you need to explain better. Explain your purpose. Explain your graph. Explain your conclusions.

Your report should contain the following sections:

*Introduction.* Discuss the goals of your kinetic analysis. Give background on copolymer kinetic analysis here, including a brief discussion of the fitting technique you plan to use. Other than to satisfy the requirements of this class, what was the aim of this experiment? (Why would anyone want to know reactivity ratios?) Give a background for the experiment by (1) discussing any theoretical issues required to understand the lab and (2) describing why the experiment might be important. If you write about something you learned from somewhere else, reference your source; however, as in journal articles, reference materials that a reader can access (e.g., other journal articles, textbooks, handbooks) rather than materials that others don't have access to (e.g., my lectures, your lab instructions). Chemical structures of your molecules, if you would like to include them, can go here as well.

It will be tempting for you to copy the Introduction of the lab instructions into the Introduction of your manuscript. They serve different purposes, and your Introduction should be more concise and general than mine anyway. So don't do this.

*Experimental.* How did you do the experiment? Your description should allow someone who has never met you or read the assignment to reproduce your work. Omit parts that don't help this imaginary person. (Example: "The aluminum pan was labeled with our names." Your imaginary experimenter will probably succeed if he labels the pan something else.) But include parts that might be important. (Example: "The sample was placed in an aluminum pan and dried in a vacuum oven for 2 days." What if the experimenter used something else other than aluminum? Heat might not transfer as well, flat shape might be important to effective curing....) As in the Experimental sections of ACS journal articles, your Experimental should be written in the third person past-tense.

Be as numerical as you can be here. Always report precisely what you did ("0.05 g (0.2 mmol) benzoyl peroxide was added") rather than what you calculated ("one molar equivalent of benzoyl peroxide was added"). In addition, be descriptive about what happened in lab. If your reaction turned blue, write that. If you departed from the written procedure, write that. It may be tempting for you to edit my lab instructions directly into your experimental; because my instructions include details that are helpful, but not really critical, you should resist this temptation and write your experimental from your lab notebook in your own words.

*Results and Discussion.* What were the results of your experiments? What did they demonstrate? Was your hypothesis proven? In this section, you should include:

- A discussion of what you did and what data you (and the class) collected.
- Any equations you need to explain your analyses.
- An example NMR spectrum with integrations, illustrating how you determined  $f$  values.
- The appropriate graph for determining  $r_1$  and  $r_2$ , (be sure to include a figure caption).
- The calculated values of  $r_1$  and  $r_2$ .
- A one-paragraph discussion of your results with a comparison to literature values. What are the sources of the error you estimated?

I have asked for the *Results and Discussion* sections together because *Macromolecules* does. This means you should present your data and discuss it at the same time. Figures should have

Figure Captions that describe the Figure and explain any peculiar notation. (If you are using Microsoft Word, you can format these automatically by selecting your Figure, and then choosing Insert -> Reference -> Caption.) If there are any sources of error that would explain why your data doesn't work the way you expected, you should include that here.

*References.* You may format your references as either footnotes or endnotes. In either case, the in-text reference number should be a superscript number. Make sure references are in ACS format: (Authors *Journal*, **Year**, *Volume*, Page) for journal articles; see the ACS Style Guide for other sources such as books.

Some resources on what the different sections of a research paper should look like are listed on the course website, on the "Links" webpage. A few tips:

- You will want to become familiar with Equation Editor (included in Microsoft Word) or another equation formatting tool (such as MathType Lite, free on the web) for your reports. It's pretty difficult to format complex equations on one line, and copying equations from my lab instructions won't format very well in your document. Help with these tools is available in Word or on the web.
- Same is true for ChemDraw or ISISDraw. The latter is a free download.
- Same is true for subscripts and superscripts. Most programs have a "Format" -> "Font" type option that will allow you to do these formats.
- Variables, and most constants, are italicized. Look at a textbook, or the ACS Style Guide, if you are unsure about whether to italicize something.
- Numbers and units are always separated by a space. ("10 mL" rather than "10mL".) So are math functions. (" $r = 1$ ", not " $r=1$ ".)
- Units of measure in ACS publications are always abbreviated. Some of these abbreviations you will feel comfortable using. (Like "mL" for milliliters.) Others you won't be. (Like "min" for minutes.) However you feel about it, abbreviate all units of measure.
- Your Experimental should be very specific. Don't make your reader calculate anything; say exactly how much of things you used, how you interpreted data, and make sure significant figures are roughly correct. ("Five milliliters" means EXACTLY five milliliters. "5 mL" means 5 +/- 1 mL, and I think you can measure more accurately than that. "5.0 mL" would probably be appropriate.)
- Do not reference my lab instructions. You will find everything that is in them in textbooks and in the papers I cite. Read, and then reference, these instead.

### **Revised Lab Report** (due in lecture, Monday, March 11)

You will have the opportunity to revise your report after we grade it the first time. You will be able to recover up to half of the maximum points for the lab report. (So if you received 70/100 on the first grade, your revision could net you 100/100; but if you received 40/100, your revision could earn you a maximum of 90/100.) If you lost points because you turned in your report late, you will *not* be able to recover those points in the rewrite. We will hand back your report in class

on March 4, with suggestions on how to improve your report. This list may not be exhaustive; if I don't circle every instance of an error, that doesn't mean I don't expect you to correct every instance of an error. You are particularly encouraged to meet with Andy Taton in his office (Smith 425) after lab to talk specifically about what you might improve in this report.

## Chem/MatS/ChEn 4223W

### Lab 2 Appendix A

#### *Processing NMR Data with ACD/NMR Processor*

It is possible to process and print NMR data in the NMR facility, but it is usually easier to work up data on your own computer, or on the computers in the Microcomputer Lab (101D Smith). In addition, when you analyze NMR data on your computer it is easier to cut and paste the spectra into your lab reports. To do this, you will need two pieces of software:

- A client program to transfer (via FTP) your NMR data from the NMR server. There are a lot of these out there. PC users can download **WinSCP** (<http://winscp.net/>) or **FileZilla** (<http://filezilla-project.org/>, also available on Microcomputer Lab workstations) for free. Mac users can also download **FileZilla** or **Cyberduck** (<http://cyberduck.ch/>) for free. These instructions will describe WinSCP (the program available in the Microcomputer Lab), but the other programs work very similarly.
- A program to Fourier transform the data and analyze the resulting spectrum. There are only a few of these. For PC users, ACDLabs has recently made **ACD/NMR Processor** free to all academic users, and this software will be described in these instructions. Browse to [http://www.acdlabs.com/resources/freeware/nmr\\_proc/index.php](http://www.acdlabs.com/resources/freeware/nmr_proc/index.php) to install the software on your own computer, or you may also use it pre-installed in the Microcomputer Lab (101D Smith). Mac users can look into **iNMR** ([www.inmr.net](http://www.inmr.net)). I think this program can do everything that you need for this course. The full version of the software costs money, but the trial version is free and appears to work the same as the full version except for occasional reminders to upgrade (which don't seem to stop you from using the software anyway). The program works much like ACD/NMR Processor, so you can use these instructions as a general guide to that software, but I won't be describing it explicitly (mainly because I've never used it).

#### **Transferring NMR data via FTP**


You will need to transfer the .fid folder from the NMR lab server to your personal computer using an FTP program such as FileZilla or WinSCP. If you use WinSCP:

- Make sure you are on-campus, or that you are using VPN tunnel software (<http://www.oit.umn.edu/vpn/>).
- Open WinSCP, and click **New** to create a new connection.
- In the next window, enter
  - Host name: echo.chem.umn.edu
  - User name: zzzlab
  - Password: *zzzlab's password*
  - File protocol: FTP




- Press **Login** to connect to the server. (You can also click **Save** to create a profile with these settings.)
- On login, the program will show the remote server's directory structure on the right, and your local computer's directory on the left. Find your `.fid` subdirectory, and click-and-drag to copy that subdirectory to your local machine. The subdirectory will contain files titled `fid.dat`, `log.dat`, `procpur.dat`, and `text.dat`; you will need them all, so just copy the directory that contains them.

## Analyzing NMR data with ACD/NMR Processor




1. Start the 1D NMR Processor program, choose to **Open an Existing File**, click **More Files...**, and browse to find your `fid.dat` file. (Make sure the File Type dropdown box lists “All Files (\*.\*)” so you can see it.)
2. The program window will show the FID:




From here you have a couple of options. If you pick the “Shortcut” button above the FID, the program will attempt to do all of the spectral analysis for you—pick peaks, find multiplets, calculate coupling constants, etc. Most of the time it will fail, either because peaks from impurities make it over-assign the spectrum, or because it attempts to classify obviously overlapping multiplets as single resonances. I think you will find that you are much more successful at it on your own.

So, instead, I recommend that you analyze your spectrum yourself. Click **Fourier Transform** to convert your FID into a 1D NMR spectrum. The program will do a few things automatically for you; it will automatically apply a phase correction, a baseline correction, and an apodization (window) function to your spectrum. In general, the software does an excellent job of all of these things, and you may never need to change the software’s choices for these. But in case you do:


- Use the zoom buttons, and particularly the Horizontal Zoom button, to focus in on an area of the spectrum that shows some structure. Click and drag to zoom in. If you want to zoom back out, click the Show Whole Spectrum button. If you have a wheel mouse, the wheel in the center can be used to control the vertical zoom independently of the horizontal zoom.






- Click **Phase** to modify the software's auto-phasing. There are two phase corrections applied, a zero-order (Ph0) and first-order (Ph1) one. The easiest way to see how this affects your spectrum is to select the **Mouse Phase** option in the command bar; in a two-button mouse, the left button controls Ph0, and the right button a combination of Ph0 and Ph1. If you don't like what you've done to the spectrum, you can always back out of the phasing menu by clicking the red X, or go all the way back the beginning with **File** → **Reopen**.
- Click **Baseline** to set baseline points; this can occasionally be useful for quantitative integration, but you probably won't ever need it.
- Click **Interactive FT** (or **WFunctions** or **Apodization** when the FID is showing) to perform resolution enhancement. Noise in the early and/or late part of the FID results in broadening of peaks in the Fourier-transformed spectrum, and so sometimes it is useful to multiply the FID by a "window" function to fix this problem.

before apodization




after apodization



In the Interactive FT window, you can perform “resolution enhancement” by switching the Window Function entries from Default to User, and then adjusting the given values for EM (exponential function) or GM (Gaussian function).

3. Back in the main spectrum window, click on **Reference** to make sure the chemical shift scale is correctly referenced to  $\text{CHCl}_3$  ( $\delta = 7.26$  ppm) and tetramethylsilane (TMS,  $\delta = 0.00$  ppm). Because the NMR instrument’s deuterium lock generally prevents the frequency from drifting, most spectra will already be extremely close to referenced. But just to make sure, move your cursor over the peak that you think to be  $\text{CHCl}_3$  or TMS, left-click, and pick the appropriate reference frequency from the list. (See next page.)



4. Next, we'll work on making the spectrum presentable. In order to put chemical shift labels on the peaks, click **PeakPicking**. If you then click **Auto** to have the computer pick peaks for you, it only does an okay job. To do a better job, zoom in on each area you want to pick peaks in, click **Peak Level**, and click your cursor at the height above which you want all peaks labeled. Or, choose **Peak by Peak** and click on peaks individually. If you don't like what you've chosen, click the **Clear** button. If you do like what you've done, click the green check mark to go back to the main spectrum window.

Your  $^1\text{H}$  peak labels should show chemical shift to the fourth decimal place. If they don't, choose **Options** → **Preferences** → **Peaks** and change the **Numerical Display** option.

5. Click **Integration** to set up integration ranges. If you pick **Auto**, the computer will not integrate intelligently—it will split up or join multiplets, and make the ranges too broad or narrow. Instead, choose **Manual**, and then click and drag regions of the spectrum you want integrated. If you click and drag a region that you want to serve as a reference—say a peak you know to correspond to a methyl  $-\text{CH}_3$ —then right after you set that integral, type an integral (3 in this case) in the **Reference** box.

If you don't provide a reference value, occasionally the program will report integrals with too few significant digits. The default setting is to give values to the hundredths place, so you will occasionally see integrals with values "0.01" or even "0.00". These values are clearly unreliable—0.01 could mean anything between 0.005 and 0.0149—but you can fix that problem by setting one of the integrals to a higher **Reference** value.

6. Most of the time, steps 1-5 will be all you need to process your spectrum. At this point, you might **File** → **Print** in order to have a hard copy, or you might **Edit** → **Copy to Clipboard** so you can paste the spectrum into your lab report.

**Lab 2 Appendix B**

*Nonlinear Least-Squares Curve Fitting with Kaleidagraph*

If you would like to try to fit your ( $f$ ,  $F$ ) data directly to the copolymerization equation, you might use Synergy Kaleidagraph. (You can also do nonlinear least-squares fitting on Microcal Origin, Mathworks MATLAB, or Wolfram Mathematica, if you have access to them and know how to use them. It is also possible to use Excel's "Solver" functionality to do nonlinear least-squares fitting, but it is complicated, and I recommend against it.) Kaleidagraph is easy to use, but you will have to go to the Microcomputer Lab (101D Smith) to use it.

To use Kaleidagraph for nonlinear least-squares fitting:

1. Open Kaleidagraph and click **File, Import, Excel** to import data from an Excel spreadsheet.
2. Graph your data by clicking **Gallery, Linear, Scatter**. Pick the appropriate columns for the x- and y-axes, and create a New Plot.
3. To fit any of your data to a custom function, click **Curve Fit, General, fit1**. (Fit1 is already defined, but we will change it. If you want to create and save functions rather than just editing the one in the program, click **Curve Fit, General, Edit General** and create new functions.)
4. Check the data box you want to fit. Then, click **Define** to create your function. The function has to have a specific form in which the independent, x-axis variable (e.g.,  $f$ ) is named "m0" and all other variables to be optimized are named "m1", "m2", etc. The equation is then followed by initial guesses for the optimized variables. So, for example, if you were graphing ( $f$ ,  $F$ ) data and you wanted to fit it to the equation

$$F_1 = \frac{r_1 f_1^2 + f_1(1 - f_1)}{r_1 f_1^2 + 2f_1(1 - f_1) + r_2(1 - f_1)^2},$$

you might type

$$\begin{aligned} & (m1*(m0)^2 + m0*(1-m0))/(m1*(m0)^2 + 2*m0*(1-m0) + m2*(1-m0)^2); \\ & m1 = 1; m2 = 1 \end{aligned}$$

in the function box. In this equation, the " $F_1 =$ " part doesn't show up, m0 stands for  $f_1$ , m1 for  $r_1$ , and m2 for  $r_2$ . You also need to input initial guesses at the end for m1 and m2; here I've guessed 1 for both. Importantly, none of your initial guesses can be set to a value of 0, or you will get an error. The program will turn your guesses into optimized, fit values for m1 and m2.

5. Click **OK**, and a fit curve should appear in your data window, along with the optimized values and other information from your fit. If you don't like your fit, you can remove it from the graph by going back into **Curve Fit, General, fit1** and deselecting the check box for the fitted data. Click **OK**, and the fit curve and fit data disappear. (This way, if you made a mistake in the equation, you can go back and re-fit onto the same data.) If you get a "Singular Coefficient Matrix" error during your fit, it is probably because your initial guesses were really far off; try new guess values.